

700575



HIGH PERFORMANCE IN CONTACTOLOGY



SL500L

DIGITAL SLIT LAMP

Range of filters



Illumination



Magnifications



Optometry





THOROUGH EXAMINATION

- High quality optical system offers wide field of view and qualitative images.
- 5 drum steps magnifications allow clear observation of finest details.
- Numerous slit adjustments and filter options provide flexibility in viewing different structures of the eye.
- Integrated yellow filter and background illumination provide optimal examination conditions.

CUTTING EDGE LED ILLUMINATION

- Brighter and more homogeneous illumination.
- Very comfortable examination for the patient with no heat emission.
- Long lifetime with substantial maintenance cost reduction.

IMAGE AND VIDEO CAPTURING (OPTION)

- Combined with the DS 550 and Anaeyes software, possibility of capturing images and video sequences.

MAGNIFICATION		SLIT	
Type	Galilean magnification changer	Slit width	Continuous from 0mm to 12mm
Standard Eyepiece	12.5x	Slit length	Continuous from 0mm to 12mm
Selection	Rotation drum 5 positions	Slit angle	0 to 180° with horizontal scanning capability
Field of view	6x (37 mm)	ILLUMINATION LAMP	
	10x (24 mm)	LED	
	16x (14 mm)	FILTERS	
	25x (8 mm)	Blue, Red-Free, Red and Yellow	
	40x (5.2 mm)	DIMENSIONS	
BASE		POWER-SUPPLY	
Longitudinal movement	113mm	Table	255(W)x 320(D) x 430(H) mm
Lateral movement	107mm	Consumption	100-120V, 220-240V AC 25 VA
Vertical movement	30mm	CE Conformity marking	
Chinrest (vertical)	76mm		
OPTIONS		DS550 Image capture device	1/1.8 CCD type Effective pixel number : 2 MP
	Polarised filters Hruby Lens Applanation tonometer		

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.

Your local distributor :